

Datasheet

MAP2 monoclonal antibody, clone MT-01

Draber P. J Histochem Cytochem. 1995
Dec;43(12):1269-84.

Catalog Number: MAB3866

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against full length native MAP2.

Clone Name: MT-01

Immunogen: Native purified bovine MAP2.

Host: Mouse

Reactivity: Bovine,Cat,Human,Mouse,Pig

Applications: ICC, IP, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody reacts with molecule Microtubule Associated Protein 2ab (MAP2ab).

Form: Liquid

Isotype: IgG1

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store at 4°C. Do not freeze.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 281294

Gene Symbol: MAP2

Gene Alias: -

References:

1. Changes of MAP2 phosphorylation during brain development. Riederer BM, Drabero E, Viklicky V,