

Datasheet

CCBE1 purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00147372-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human CCBE1 protein.

Immunogen: CCBE1 (AAH46645, 1 a.a. ~ 215 a.a)
full-length human protein.

Sequence:

MVKAGTCCATCKEFYQMKQTVLQLKQKIALLPNNAAD
LGKYITGDKVLASNTYLPGLPGGQGPSPGPKG
SPGFPGMPGPPGQPGPRGSMGPMGPSPDLSHIKQG
RRGPVGPPGAPGRDGSKGERGAPGPRGSPGPPGSF
DFLLMLADIRNDITELQEKVFGHRTSSAEFPLPQEF
PSYPEAMD LGSGDDHPRRTETRD LRAPRDFYP

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 147372

Gene Symbol: CCBE1

Gene Alias: FLJ30681, MGC50861