

Datasheet

ZNF548 MaxPab mouse polyclonal antibody (B01)

Gene Symbol: ZNF548

Gene Alias: FLJ32932

Catalog Number: H00147694-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF548 protein.

Immunogen: ZNF548 (NP_690873, 1 a.a. ~ 533 a.a) full-length human protein.

Sequence:

MNLTEGRVVFEDVAIYFSQEEWGHLEAQRLLYRDVM
LENLALLSSLGSHWGAEEDESPSQGGFSVGVSEVTTS
KPCLSSQKVHPSETCGPPLKDILCLVEHNGIHPEQHIYI
CEAELFQHPKQKQIGENLSRGDDWIPSGKNHRVHMA
EEIFTCMEGWKDLPATSCLLQHQPQSEWKPYRDTE
DREAFQTGQNDYKCSECGKTFTCSYSFVEHQIHTGE
RSYECNKGKFFKYSANFMKHQTVHTSERTYECREC
GKSFMYNYRLMRHKRVHTGERPYECNTCGKFFRYSS
TFVRHQRVHTGERPYECRECGKFFMDSSTLIKQHRVH
TGERPYKCNDGKFFRYISTLIRHQRIHTGERPYECV
CGELFRYNSSLVKHWRNHTGERPYKCSECGKSFYH
CRLIRHQRVHTGERPYECSECGKFFRYNSNLIKHW
HTGERPYECRECGKAFSHKHILVEHQIHSGERPYEC
SECQKAFIRKSHLVHHQKIHSEERLVCSMNVGNSLAKT
PTSLNIRDFTMEKVYH

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 147694