

Datasheet

CENTB2 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00023527-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human CENTB2 protein.

Immunogen: CENTB2 (NP_036419.2, 1 a.a. ~ 778 a.a) full-length human protein.

Sequence:

MKMTVDFEELCKDSPRFRAALEEVEGDVAEELKLDK
LVKLCIAMIDTGKAFVCVANKQFMNGIRDLAQYSSNDV
VETSLTKFSDSLQEMINFHTILFDQTQRSIKAQLQNFVK
EDLRKFDAKKQFEKVSEEKENALVKNAQVQRNKQH
EVEEATNILTATRKCFRHALDYVLQINVLSKRRSEILK
SMLSFMAYHLAFFHQGYDLFSELGPYMKDLGAQLDRL
VVDAAKEKREMEQKHSTIQKDFSSDDSKLEYNVDAA
NGIVMEGYLFKRASNAFKTWNRRWFISIQNNQLVYQKK
FKDNPTVVVEDLRLCTVKHCEDIERRFCFEVVSPTKSC
MLQADSEKLRQAWIKAVQTSIATAYREKGESEKLDK
KSSPSTGSLDSGNESKEKLLKGESALQRVQCIPGNAS
CCDCGLADPRWASINLGITLCIECSGIHRS LGVHFSKV
RSLTLDTWEPELLKLMGELGNDVINRVYEANVEKMGIK
KPQPGQRQEKEAYIRAKYVERKFVDKYSISLSPPEQQ
KKFVSKSSEEKRLSISKFGPGDQVRASQSSVRSNDS
EIQSSDDGRESLPSTVSANSLEYEPGERQDSSMFLD
SKHLNPGQLYRASYEKNLPKMAEALAHGADVNWAN
SEENKATPLIQAVLGGSLVTCEFLQNGANVNQRDVQ
GRGPLHHATVLGHTGQVCLFLKRGANQHATDEEGKD
PLSIAVEAANADIVTLLRLARMNEEMRESEGLYGQPGD
ETYQDIFRDFSQMASNNPEKLNRFQQDSQKF

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 23527

Gene Symbol: ACAP2

Gene Alias: CENTB2, CNT-B2, KIAA0041