

## Datasheet

### KLF15 monoclonal antibody, clone 2G8

**Catalog Number:** MAB10713

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against partial recombinant KLF15.

**Clone Name:** 2G8

**Immunogen:** Recombinant protein corresponding to human KLF15.

**Host:** Mouse

**Theoretical MW (kDa):** 44

**Reactivity:** Human

**Applications:** ELISA, WB-Ce

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Isotype:** IgG1

**Recommend Usage:** ELISA (1:10000)

Western Blot (1:500-1:2000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In ascites (0.03% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 28999

**Gene Symbol:** KLF15

**Gene Alias:** DKFZp779M1320, KKLF