

Datasheet

KLHL1 monoclonal antibody, clone 3A8

Catalog Number: MAB10716

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant KLHL1.

Clone Name: 3A8

Immunogen: Recombinant protein corresponding to human KLHL1.

Host: Mouse

Theoretical MW (kDa): 82

Reactivity: Human

Applications: ELISA, IF, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Isotype: IgG1

Recommend Usage: ELISA (1:10000)
Western Blot (1:500-1:2000)
Immunofluorescence (1:200-1:1000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In ascites (0.03% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 57626

Gene Symbol: KLHL1

Gene Alias: FLJ30047, KIAA1490, MRP2