

Datasheet

KLHL11 monoclonal antibody, clone 1B9C1

Gene Alias: FLJ10572**Catalog Number:** MAB10718**Regulation Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody raised against partial recombinant KLHL11.**Clone Name:** 1B9C1**Immunogen:** Recombinant protein corresponding to human KLHL11.**Host:** Mouse**Theoretical MW (kDa):** 80**Reactivity:** Human**Applications:** ELISA, IF, IHC-P, WB-Ce
(See our web site product page for detailed applications information)**Protocols:** See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Form:** Liquid**Isotype:** IgG1**Recommend Usage:** ELISA (1:10000)
Western Blot (1:500-1:2000)
Immunohistochemistry (1:200-1:1000)
Immunofluorescence (1:200-1:1000)
The optimal working dilution should be determined by the end user.**Storage Buffer:** In ascites (0.03% sodium azide)**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.**Entrez GeneID:** 55175**Gene Symbol:** KLHL11