

Datasheet

KLHL12 monoclonal antibody, clone 2G2

Catalog Number: MAB10719

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant KLHL12.

Clone Name: 2G2

Immunogen: Recombinant protein corresponding to human KLHL12.

Host: Mouse

Theoretical MW (kDa): 63

Reactivity: Human

Applications: ELISA, WB-Ce
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Isotype: IgG1

Recommend Usage: ELISA (1:10000)
Western Blot (1:500-1:2000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In ascites (0.03% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 59349

Gene Symbol: KLHL12

Gene Alias: C3IP1, DKIR, FLJ27152