

Datasheet

FGF21 purified MaxPab rabbit polyclonal antibody (D01P)

and invasion. The function of this growth factor has not yet been determined. [provided by RefSeq]

Catalog Number: H00026291-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human FGF21 protein.

Immunogen: FGF21 (AAH18404.1, 1 a.a. ~ 209 a.a) full-length human protein.

Sequence:

MDSDETGFEHSGLWVSVLAGLLLGACQAHPIPDSSPL
LQFGGQVRQRYLYTDDAQQTEAHLEIREDDGTGGAA
DQSPELLQLKALPGVILGVKTSRFLCQRPDGALY
GSLHFDPEACSFRELLLEDGYNVYQSEAHGLPLHLP
NKSPHRDPAPRGPARFLPLPGLPPAPPEPPGILAPQP
PDVGSSDPLSMVGPSQGRSPSYAS

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 26291

Gene Symbol: FGF21

Gene Alias: -

Gene Summary: The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities and are involved in a variety of biological processes including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth