

Datasheet

MYCBP MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00026292-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MYCBP protein.

Immunogen: MYCBP (AAH08686, 1 a.a. ~ 103 a.a)
full-length human protein.

Sequence:

MAHYKAADSKREQFRYLEKSGVLDLTCLKVLVALYEE
PEKPNALDFLKHHLGAATPENPEIELLRLELAEMKEK
YEAIVEENKKLKAKLAQYEPPEEKRAE

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 26292

Gene Symbol: MYCBP

Gene Alias: AMY-1

Gene Summary: The MYCBP gene encodes a protein that binds to the N-terminal region of MYC (MIM 190080) and stimulates the activation of E box-dependent transcription by MYC.[supplied by OMIM]