

## Datasheet

### TAS2R49 polyclonal antibody

**Catalog Number:** PAB25086

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of TAS2R49.

**Immunogen:** A synthetic peptide corresponding to TAS2R49.

**Host:** Rabbit

**Theoretical MW (kDa):** 38

**Reactivity:** Human

**Applications:** IF, WB-Ce

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** TAS2R49 polyclonal antibody detects endogenous levels of TAS2R49 protein.

**Form:** Liquid

**Purification:** Affinity purification

**Concentration:** 1 mg/mL

**Recommend Usage:** Western Blot (1:500-1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.2 (0.05% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 259295

**Gene Symbol:** TAS2R49

**Gene Alias:** MGC119118, MGC119119, MGC119120, MGC119121, MGC134637, T2R49, T2R56