

Datasheet

DPYSL4 monoclonal antibody (M01), clone 1F5

Catalog Number: H00010570-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant DPYSL4.

Clone Name: 1F5

Immunogen: DPYSL4 (NP_006417, 461 a.a. ~ 559 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

VTPGAGRFVPRKTFPDFVYKRIKARNRLAEIHGVPRGL
YDGPVHEVMVPAKPGSGAPARASCPGKISVPPVRNLH
QSGFSLSGSQADDHIARRTAQKIM

Host: Mouse

Reactivity: Human

Applications: ELISA, IP, S-ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10570

Gene Symbol: DPYSL4

Gene Alias: CRMP3, DRP-4, ULIP4