

Datasheet

TINF2 monoclonal antibody (M02), clone 3G11

Catalog Number: H00026277-M02

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant TINF2.

Clone Name: 3G11

Immunogen: TINF2 (NP_036593, 256 a.a. ~ 354 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

RHFNLA PLGRRRVQSQWASTRGGHKERPTVMLFPFR
NLGSPTQVISNPESKEEHAIYTADLAMGTRAPSNQKYK
GPYQTLGGRALKENPVDLPATEQKE

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 26277

Gene Symbol: TINF2

Gene Alias: TIN2, TIN2L